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Addition-vulcanising vinyl-polysiloxane in several final hardnesses for all dental duplications. Easy and precise to mix in a ratio of 1:1.

Product description

Dublisil®-products are universal duplication materials for all kinds of duplications. Due to the practical way of mixing and the very thin-flowing initial consistency absolutely price duplication models can be produced. Dublisil®-forms are nearly shrinkless and have a high form consistency as well as a long shelf life due to the used platinum hardener system. Negative form made of Dublisil® can easily and without losing precision be poured out for several times The surface of the duplication model is smooth and hard without after-treatment.

Preparatory measures

Almost every material e. g. plaster models, model resins and metals can be duplicated with Dublisil®. Important: The master model must not be watered but should be dry and free of any insulating agents etc. For an economical consumption of the duplication material we recommend to make use of the Dublisil®-flask (D3805) with a variably adjustable interior as well the foam space maintainer (D38076) for the blocking out of the base parts which are unimportant for the duplication.

Dosage

Dublisil® is mixed in the ratio 1:1 of the components A (catalyst / white) and B (base / coloured).

Processing

The two single components can manually be mixed in a mixing cup or automatically by means of a dosage and mixing unit. The most rational and clean method is the mixing in the Dosper evo. At the touch of the button Dublisil® is fully automatic and losslessly measured out, mixed and produced without air bubbles. In case that the material is manually mixed the fill in the material in a preferably thin stream into the flask. When using the Dosper evo this stream is defined by the unit. The distance to the flask should be approx. 20 cm.

Further processing

The best results regarding precision and homogeneity are achieved by the so-called pressure duplication technique. For this the just filled flask is put into a pressure pot (e. g. Polymax). In case that the duplication form is produced under pressure it is important that also the model is produced under the same pressure. After the vulcanisation has been completed the master model can directly be deflasked.

General aspects

Dublisil® is unproblematic in physiological and toxicological respects. A chemical and thermal disinfection is possible. Wear protective clothing! Silicone impurities on clothes cannot be cleaned – also not chemically. Allergic reaction to product components cannot be excluded in case of sensitive persons. The vulcanisation can be disturbed by latex gloves, hand cremes, cleaning agents etc. The product is exclusively designed for the a. m. field of application and is only to be used by qualified and instructed persons.

The general information of this product is based on detailed application practice and intensive research work. Due to further development we reserve the right of modification.

Keep away from children! Only for the a. m. field of application!

Storage conditions
At room temperature, i. e. 18–28 °C. **Durability:** see date of expiry. Dublisil® is not to be used when expiry date has passed.

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This product was developed for dental use and has to be applied according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according to suitability and possible use in self dependence before using it subject to its purpose, especially in case the purpose is not mentioned in the instructions.

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